



"Our restorations aren't just **Made In The USA**. The **materials & equipment** we purchase are of the highest quality from the **most reputable brands on the planet!**"



Digital Impression Guide.

We are working with all IOS platforms, so our customers can take advantage of TrueFit's speed, accuracy and Ai enhancement features.



3Shape Trios

Use the "add connections" feature to send request to ios@truefit.digital

We will get the request within 5 minutes. Once approved the doctor will be able to send scans.



iTero

Search for "TrueFit Digital Lab" or by company ID (573253) through "add preferred lab". If either option is unavailable, please call iTero directly to add TrueFit Digital Labs as a preferred lab.

When the doctor scans anycase, the file goes to Align Tech and iTeroteam will evaluate the scan. If the scan needs to be taken again, iTerowill contact the doctor. If the scan is approved, it will be sent to the lab within 8 hours.



Shining 3D

Use "Connections" to access" + NEW" tab.

Search Truefit Digital Laboratory" via Zip Code: 83709



Allied Star

*currently being updated



Dentsply Sirona

*currently being updated



Dexis

*currently being updated



E4D

*currently being updated



Ivoclar

*currently being updated



Medit Link

*currently being updated